

**MOD'ALLI NURSING CART WITH DRAWERS**

**Description**

DIMENSIONS: L70,8 X H151,1 X D50,5 CM -MOD'ALLI LOW CABINET F600 15G WITH RETRACTABLE WRITING SHELF - INCLUDING: 1 WHITE NARCOTICS DRAWER WITH DOOR AND MECHANICAL LOCK HEIGHT 075 (4 DRAWER DIVIDERS 1/4) - 1 WHITE DRAWER HEIGHT 155 (6 DRAWER DIVIDERS 1/6) - PANEL MEDICAL LAPTOP 22 INCHES WATERPROOF IP65 - BATTERY 85 WATTS - KEYBOARD + MOUSE IP68

