



**BEDS-FIXED CLINIQUE ON FOOT PADS - FLAT BED LOW**

**Description**



- 1 - Comfort:** The bed is ergonomic and anatomical.
- 2 - EasyClean:** Its tubular structure, easy to clean, avoids any risk of retention. It saves time for decontamination.

**WEIGHT**

Weight bare structure 43 kg

**DIMENSIONS**

Height Sleeping surface from the grow : 45 cm

**SAFETY CHARGE**

Maximum patient weight 135 kg  
SWL 170 kg



## WEIGHT

Weight bare structure 43 kg

## DIMENSIONS

Length bare structure 208 cm  
 Width bare structure 90,2 cm  
 Length structure + H&F Boards + siderails 208 cm  
 Width structure + H&F Boards + siderails 99 cm  
 Length sleeping surface 200 cm  
 Width sleeping surface 90 cm  
 Height Sleeping surface from the grow : 45 cm

## TYPE

Variable height type F-Type  
 Type of engine Any

## WEIGHING SYSTEM

Weighed (Yes/No) No

## MATERIAL & COATING

Tray Structure Gray Epoxy Steel  
 Bedbase / sleeping surface Welded metal slats

## BASE

Type of footing Tubular on pads

## SAFETY CHARGE

Maximum patient weight 135 kg  
 SWL 170 kg

## CONFIGURATION

Configuration Siderails - H&F boards  
 Number of lifting pole support 2  
 Number of IV pole support 2

## REGULATION

Medical device class Class I  
 CE marking CE

## FIRE CLASSIFICATION

Melaminated panel M3

## OPTIONS

BED LOW ON 4 CASTORS DIAMETER 125 MM  
INCLUDING 2 WITH BRAKES  
**Réf. OP-4221-ROUES**



## OPTIONS

ALTERNATIVE OPTION ON PACK - BED DELIVERY  
ON VERTICAL TRANSPORT TROLLEY  
**Réf. OPP-KIT-LIV-VER**



## OPTIONS

MANUAL BACKREST ON RACKS (1P)  
**Réf. OP-4221-1P**



BACKREST ASSISTED BY JACK GAS (1PA)  
**Réf. OP-4221-1PA**



MANUAL BACKREST & LEG REST ON RACKS (2P)  
**Réf. OP-4221-2P**



BACKREST ASSISTED BY JACK GAS & MANUAL  
LEG REST (2PA)  
**Réf. OP-4221-2PA**



## ACCESSORIES

ALTERNATIVE OPTION ON PACK - BED DELIVERY  
ON VERTICAL TRANSPORT TROLLEY  
**Réf. ACP-4997-00013**



## ACCESSORIES

FOLDING BARRIERS WITH HORIZONTAL INTERLOCKING (CAM CLAMPING) - PER PAIR  
**Réf. AC-4481-00002**



SET OF FOLDABLE SIDE RAILS WITH  
HORIZONTAL FITTING (CLAMPS TIGHTENING)  
(THE PAIR)  
**Réf. AC-4481-00003**



SET OF TWINGUARD FOLDABLE STEELIC SIDE  
RAILS WITH HORIZONTAL FITTING  
**Réf. AC-4481-00004**



GRIPPING HANDLE - EASE GETTING OUT OF BED - CLAMPS TIGHTENING (UNIT)  
**Réf. AC-4996-00067**



SET OF PADDED PROTECTIONS FOR CLASSIC  
SIDE RAILS (THE PAIR)  
**Réf. AC-LX-EVO07-2440P**



## ACCESSORIES

LIFTING POLE - EPOXY FINISH  
**Réf. AC-4300-00001**



LIFTING POLE WITH REEL STRAP - EPOXY FINISH  
**Réf. AC-4300-00002**



OPTION 4 HOOKS ON IV POLE ON WHEELS  
**Réf. OP-4996-4CRO**



## ACCESSORIES

ADDITIONAL SIDE RAILS TO PLACE AT THE  
FOOT OF THE BED - (THE PAIR)  
**Réf. AC-4996-00090**



**ACCESSORIES**

IV POLE WITH 2 HOOKS  
**Réf. AC-4996-00006**



IV POLE (4 HOOKS)  
**Réf. AC-4996-00031**



TELESCOPIC IV POLE WITH 2 HOOKS  
**Réf. AC-4996-00032**



TELESCOPIC IV POLE (4 HOOKS)  
**Réf. AC-4996-00033**



**ACCESSORIES**

MOBILE IV POLE WITH 2 HOOKS  
**Réf. AC-4996-00034**



**ACCESSORIES**

BED LOW EXTENSION 25 CM REMOVABLE & NON ADJUSTABLE FOR BED 90 CM  
**Réf. AC-4996-00028**



BED LOW EXTENSION ADJUSTABLE 25 CM  
**Réf. AC-4996-00029**



EXTENTION FOR WAFFLE FOAM MATTRESS + COVER - 20 X 90 X 14 CM  
**Réf. AC-LX-M315-Y4**



MATTRESS + COVER - VISCOELASTIQUE - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION  
**Réf. AC-LX-M325-Y4**



**ACCESSORIES**

MATTRESS + COVER - COMFORT - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION  
**Réf. AC-4000-00036**



COVER FOR MATTRESS - WIDTH 90 CM - TRUNCATED ANGLES  
**Réf. AC-LX-H007-A-Y4**



COVER FOR MATTRESS - WIDTH 90 CM - RIGHT ANGLES  
**Réf. AC-LX-H007-Y4**



MATTRESS + COVER - COMFORT MULTI-ARTICULATION - SLEEPING 200 X 90 CM - TRUNCATED ANGLES  
**Réf. AC-LX-M207-A-Y4**



MATTRESS + COVER - MULTI-ARTICULATION COMFORT - LYING SURFACE 200 X 90 CM - RIGHT ANGLES  
**Réf. AC-LX-M207-Y4**



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION  
**Réf. AC-LX-M305-Y4**



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 200 X 90 CM - TRUNCATED ANGLES  
**Réf. AC-LX-M407-A-Y4**



MATTRESS + COVER - MULTI-JOINTS WAFFLE - LYING SURFACE 200 X 90 CM - RIGHT ANGLES  
**Réf. AC-LX-M407-Y4**



MATTRESS + COVER - VISCOELASTIC - SLEEPING 200 X 90 CM - TRUNCATED ANGLES  
**Réf. AC-LX-M507-A-Y4**



MATTRESS + COVER - VISCOELASTIC - LYING SURFACE 200 X 90 CM - RIGHT ANGLES  
**Réf. AC-LX-M507-Y4**



**ACCESSORIES**

SHOCK ABSORPTION MAT - PROTECTION IN CASE OF FALL  
**Réf. AC-LX-MS200**



PATIENT PRESENCE AND/OR POSITION SENSOR IN THE BED  
**Réf. AC-LX-SANTINELLP**



**ACCESSORIES**

UNIVERSAL LOWIN HOLDER (ADAPTATION TO HEAD AND FOOT BOARDS)

**Réf. AC-4996-00003**



URINAL POCKET HOLDER

**Réf. AC-4996-00018**



URINAL HOLDER

**Réf. AC-4996-00045**



URINAL BAG HOLDER WITH WALL-MOUNTED SUPPORT

**Réf. AC-4996-00046**



LOWSIN HOLDER

**Réf. AC-4996-00120**

