

BED PACK Y-TYPE Y3 FRONT-EXIT MOVE WIDTH 90 CM - CLASSIC SIDE RAILS AND SPECIFIC MATTRESS



Description

1 - EasyNursing: The H-frame has been specifically designed to be lighter, with no retention zone. With its full-width sleeping surface, patient handling and comfort are optimised. We must work with you to reduce MSD.

2 - Comfort: The Front Exit function offers the patient a real sitting position, humanizing during tours, suitable for meals and helping to get up.

3 - Ultra low: The Y-Type Ymagyn 4 bed offers a very low position to secure the patient in the event of a fall during the night. Restraint can thus be avoided.

4- Full Accessorization: Its EasyClip system offers a wide choice of accessories in order to offer a bed adapted to each pathology.

WEIGHT

Weight bare structure	90 kg
-----------------------	-------

DIMENSIONS

Length sleeping surface	200 cm
Width sleeping surface	90 cm

TYPE

Variable height type	Y-TYPE
----------------------	--------

MATERIAL & COATING

Bedbase / sleeping surface	Compact HPL
----------------------------	-------------

SAFETY CHARGE

Maximum patient weight	195 kg
SWL	260 kg

WEIGHT	
Weight bare structure	90 kg
DIMENSIONS	
Length bare structure	210 cm
Width bare structure	89,7 cm
Length structure + H&F Boards + siderails	215 cm
Width structure + H&F Boards + siderails	100 cm
Length sleeping surface	200 cm
Width sleeping surface	90 cm
Variable height minimum	20 cm
Variable height maximum	80 cm
Automatic medium stop position	40 cm
Automatic security stop position	24 cm
TYPE	
Variable height type	Y-TYPE
Type of engine	Electric
Type of brake	Centralized brake in foot only
WEIGHING SYSTEM	
Weighed (Yes/No)	No
MATERIAL & COATING	
Tray Structure	White Epoxy Steel
Base Structure	Black Epoxy Steel
Variable height	Black Epoxy Steel
Bedbase / sleeping surface	Compact HPL
Wheels	Thermoplastic rubber, grey non-marking
BASE	
Type of footing	Castors
Wheels diameter	75 mm

INCLINATIONS	
Backrest	74 °
Knee brake section	38 °
Knee brakes	130 °
Translation	7 cm
Trendelenburg	17 °
Reverse trendelenburg	15 °
PROTECTION SIGN	
Bed	IPX4
Engines	IP66
SAFETY CHARGE	
Maximum patient weight	195 kg
SWL	260 kg
CONFIGURATION	
Configuration	Siderails - H&F boards - Mattress
REGULATION	
Medical device class	Class I
CE marking	CE
European regulation (UE 2017/745)	Yes
FIRE CLASSIFICATION	
Motorization	V0
Compact HPL	M1
Laminated	M1
Melaminated panel	M3

OPTIONS

EASYTOUCH CONTROL SATELLITE WITH PRE-PROGRAMMING OF THE MAIN FUNCTIONS INSTEAD OF STANDARD REMOTE-CONTROL
Réf. OP-L-465-ETTILH



EASYTOUCH CONTROL SATELLITE WITH PRE-PROGRAMMING OF THE MAIN FUNCTIONS WITH REMOTE-CONTROL
Réf. OP-L-465-ETTILHTEL



BACKUP BATTERY
Réf. OP-L-470-B-Y4



INTERNAL EQUALISER LOCATED ON HEAD PART SECTION
Réf. OP-L-486-EQUI-Y



NIGHT LIGHT PLACED UNDER THE SLEEPING PLATFORM
Réf. OP-L-550-LU



FUNCTION CONDEMNATION BOX WITH REMOTE-CONTROL OUTLET PLUG & PLAY FOR BED FRONT EXIT
Réf. OP-L-CONDAMNAT-FE



OPTIONS

BILATERAL RELEASE HANDLE OF BACKREST - FOR BED FRONT EXIT
Réf. OP-L-314



BILATERAL RELEASE HANDLE OF BACKREST - FOR BED FRONT EXIT
Réf. OP-L-314-2



OPTIONS

ALTERNATIVE OPTION ON PACK - MOBILITY OF THE BED Y4 EMOVE
Réf. OP-L8406-EMOVE



OPTIONS

ALTERNATIVE OPTION ON PACK - SET OF H&F BOARDS CINTEA - STEEL HANGER - HORIZONTAL FITTING
Réf. OP-L-750-CINTEA-Y4



2 ANNULAR PROTECTIVE BUMPERS FIXED ON THE HEAD & FOOT PANELS
Réf. OP-L-860-BU-ANN



OPTION HANDRAIL AT THE HEAD FOR SET OF ROUND H&F BOARDS VITTORIA
Réf. OP-LPAN-MC



ALTERNATIVE OPTION ON PACK - SET OF ONE-PIECE ROUND H&F BOARDS ORELINE Y4
Réf. OP-LPAN01-OR



ALTERNATIVE OPTION ON PACK - SET OF SQUARED WOODEN H&F BOARDS VITTORIA
Réf. OP-LPAN12-VC



ALTERNATIVE OPTION ON PACK - SET OF UPHOLSTERED ROUND H&F BOARDS CAPRI
Réf. OP-LPAN30-CR



ALTERNATIVE OPTION ON PACK - ROUND H&F BOARDS VITTORIA ***LESS VALUE***
Réf. OP-MV-DOSSERETS



ALTERNATIVE OPTION ON PACK - H&F BOARDS PPS - TWO POSITIONS BOARD AND TRAY
Réf. OP-PPS-Y4



OPTIONS

CHOCOLATE FINISH OF THE 2 TOP BARS FOR SET OF CLASSIC FOLDABLE SIDE RAILS 3 BARS - FOR BEDS FRONT EXIT
Réf. OP-L-800-HB-CHO



NATURAL FINISH OF THE 2 TOP BARS FOR SET OF CLASSIC FOLDABLE SIDE RAILS 3 BARS - FOR BEDS FRONT EXIT
Réf. OP-L-800-HB-NAT



SET OF CLASSIC FOLDABLE SIDE RAILS - 3 BARS FOR BED Y3 FRONT EXIT
Réf. OP-L-813



SET OF CLASSIC FOLDABLE SIDE RAILS - 4 BARS - ADAPTED TO GROWN CHILDREN
Réf. OP-L-814



OPTIONS

4 ACCESSORY RAILS FOR TOP BED FRAME (2 ON EACH SIDE OF THE BED)
Réf. OP-L-860-RAILCAH



ACCESSORIES

FLEXIBLE HOLDER FOR WIRED REMOTE-CONTROL
Réf. AC-LX-ST2H-Y4



ACCESSORIES

SET OF PADDED PROTECTIONS FOR CLASSIC SIDE RAILS (THE PAIR)
Réf. AC-LX-EVO07-2440P



SET OF ADDITIONNAL TOP BAR FOR CLASSIC FOLDABLE SIDE RAILS (THE PAIR)
Réf. AC-LX-REHAUSSE-B



ACCESSORIES

SET OF SMARTGUARD 60/40 TOTAL CLOSURE - WOODEN FINISH - WITHOUT H&F BOARDS - FOR BED FRONT EXIT
Réf. AC-L-840-SG-B-FE



SET OF PADDED PROTECTIONS FOR SMARTGUARD SIDE RAILS
Réf. AC-LX-PROTEC-BAR-SG



SET OF ADDITIONNAL TOP BAR FOR SMARTGUARD
Réf. AC-LX-REHAUSSE-SG



ACCESSORIES

FIXED LIFTING POLE WITH AUTOMATIC REEL
Réf. AC-LX-P1010-ENR



LIFTING POLE SWIVELLED BY A PULL WITH AUTOMATIC REEL
Réf. AC-LX-P1040-ENR



SUPPORT 4 PLASTIC HOOKS FOR LIFTING POLE
Réf. AC-LX-PCRO



ACCESSORIES

HANDLE TO EASE GETTING OUT OF BED (SOLD PER UNIT) - INCOMPATIBLE WITH KIT XL 110 CM (UNIT)
Réf. AC-LX-ASL-Y4



HANDLES TO EASE GETTING OUT OF BED FOR BED FRONT-EXIT (THE PAIR) - WITHOUT CUSHIONS
Réf. AC-LX-CALFEH



LATERAL PROTECTIVE CUSHIONS FOR ACCESSORY LX-CALFEH FOR BED FRONT-EXIT (THE PAIR)
Réf. AC-LX-COUSCALFEH



HANDLES TO EASE GETTING OUT OF BED FOR BED FRONT-EXIT (THE PAIR)
Réf. AC-LX-PALFEH



ACCESSORIES

IV POLE WITH FIXED HEIGHT WITH 4 HOOKS
Réf. AC-LX-TPS-4C-Y4



IV POLE WITH HEIGHT ADJUSTABLE BY SCREW SYSTEM WITH 4 HOOKS
Réf. AC-LX-TPSR



ACCESSORIES

COMFORTGUARD (UNIT)
Réf. AC-LX-ACG



ACCESSORIES

COVER FOR MATTRESS - WIDTH 90 CM - FOR FRONT-EXIT BED
Réf. AC-HOUSSEFEH



MATTRESS + COVER - ANTI-BEDSORE - SLEEPING 200 X 90 CM - FOR FRONT-EXIT BED
Réf. AC-LX-MATFEH



ACCESSORIES

SHOCK ABSORPTION MAT - PROTECTION IN CASE OF FALL
Réf. AC-LX-MS200



PATIENT PRESENCE AND/OR POSITION SENSOR IN THE BED
Réf. AC-LX-SANTINELLP



ACCESSORIES

BUMPER ON BED LOW FOR BED Y4 (HARMONY, MOVE, EMOVE)
Réf. AC-LX-BU-YHE



BUMPER ON TOP BED FRAME Y4
Réf. AC-LX-BU-YP-Y4



ACCESSORIES

MULTI-PURPOSE LOWKET HOLDER
Réf. AC-LX-APMUH



URINAL HOLDER ON SIDE RAILS
Réf. AC-LX-PU3



URINAL POCKET HOLDER
Réf. AC-LX-SPU1H



ACCESSORIES

TRANSPORT JIG TROLLEY FOR YMAGYN 3 & YMAGYN 4
Réf. AC-LX-CT-YY3



ACCESSORIES

AFTER SALES SERVICE KIT FOR YMAGYN BED VERSION FRONTXIT
Réf. SE-KT-ELECT-YMA-FE

