



OVERBED TABLE - RECLINING LARGE TOP TRAY + LATERAL PLASTIC COMPARTMENT

Description

1 - EasyClean : Its simple-to-clean tubular structure avoids any risk of retention. The side tray is composed of a removable thermoformed tray for maximum hygiene.

2 - Confort : Its variable height is easily adjusted by a simple action on the handle. The table is then adjusted to the desired position. Its reclining large tray equipped with a rack to prevent objects from falling.

3 - Security : Son Its wide footings with its double castors wheels offer a better stability. The automatic winding of the tray avoids the risk of accidental entrapment for the patient.

WEIGHT

Weight	18,5 kg
--------	---------

DIMENSIONS

Width	96,2 cm
Depth	45,5 cm

TYPE

Variable height type	Gaz spring
----------------------	------------

WEIGHT		FIRE CLASSIFICATION	
Weight	18,5 kg	ABS Plastic	M2
DIMENSIONS		ADDITIONAL INFORMATION	
Width	96,2 cm	Comment	1 book stand
Depth	45,5 cm		
Width large tray	75 cm		
Depth large tray	45 cm		
Width small tray	20 cm		
Depth small tray	45 cm		
Variable height minimum	75 cm		
Variable height maximum	110 cm		
TYPE			
Variable height type	Gaz spring		
MATERIAL & COATING			
large tray	Melamine faced particle board - Thickness 18 mm		
small tray	Removable ABS tray		
Column Structure	Extruded aluminum		
Tray Structure	Epoxy painted steel		
Base Structure	Epoxy painted steel		
Edges	Assoted ABS		
Wheels	Wheel centre made of Polypropylene - Thermoplastic rubber, grey non-marking		
BASE			
Type of footing	"C" base - 4 swivel castors		
Wheels diameter	50 mm		
INCLINATIONS			
Main tray	6-position rack: 0 to 42°		
SAFETY CHARGE			
Maximum admissible load	25 kg		
FIRE CLASSIFICATION			
Laminated	M1		