



VARIABLE HEIGHT OVERBED TABLE - JACK GAS ASSISTED

Description

- 1 - EasyClean :** Its simple-to-clean tubular structure avoids any risk of retention. The side tray is composed of a removable thermoformed tray for maximum hygiene.
- 2 - Confort :** Its variable height is easily adjusted by a simple action on the handle. The table is then adjusted to the desired position. Its reclining large tray equipped with a rack to prevent objects from falling.
- 3 - Security :** Son Its wide footings with its double castors wheels offer a better stability. The automatic winding of the tray avoids the risk of accidental entrapment for the patient.

WEIGHT

Weight	13 kg
--------	-------

DIMENSIONS

Width	81 cm
Depth	44 cm

SAFETY CHARGE

Maximum admissible load	5 kg
-------------------------	------

WEIGHT

Weight 13 kg

DIMENSIONS

Width 81 cm
 Depth 44 cm
 Width large tray 55,5 cm
 Depth large tray 38,5 cm
 Width small tray 25 cm
 Depth small tray 38,5 cm
 Variable height minimum 79 cm
 Variable height maximum 104,5 cm

TYPE

Variable height type Vérin à gaz

MATERIAL & COATING

Large tray Melamine faced particle board - Hetre Naturel
 Small tray Melamine faced particle board - Hetre Naturel - With fall protection gallery
 Base Structure Powder coated steel Aluminum gray
 Wheels Wheel centre made of Polypropylene - Thermoplastic rubber, grey non-marking

BASE

Type of footing "C" base - 4 swivel castors with 2 brakes
 Wheels diameter 50 mm

INCLINATIONS

Main tray Rack: 0 to 70°

SAFETY CHARGE

Maximum admissible load 5 kg

ADDITIONAL INFORMATION

Comment 1 book stand