



BED PACK Y-TYPE TIVOLI 4 - 3PEE 4 ACTUATORS MOVE - WIDTH 90 CM - H&F BOARDS VERONE - SIDERAILS CLASSIC

Description

- 1 - EasyNursing:** The H-frame has been specifically designed to be lighter, with no retention zone. With its full-width sleeping surface, patient handling and comfort are optimised. We must work with you to reduce MSD.
- 2 - Ultra low:** The Y-Type Tivoli 4 bed offers a very low position to secure the patient in the event of a fall during the night. Restraint can thus be avoided.
- 3 - Full Accessorization:** Its EasyClip system offers a wide choice of accessories in order to offer a bed adapted to each pathology.
- 4 - EasyClean:** Cleaning is made easier thanks to its structure with zero visible cables and zero retention zones.

WEIGHT

Weight bare structure	90 kg
-----------------------	-------

DIMENSIONS

Length sleeping surface	197,7 cm
Width sleeping surface	85,8 cm

TYPE

Variable height type	Y-Type
----------------------	--------

MATERIAL & COATING

Bedbase / sleeping surface	Compact HPL
----------------------------	-------------

SAFETY CHARGE

Maximum patient weight	195 kg
SWL	260 kg



WEIGHT	
Weight bare structure	90 kg
DIMENSIONS	
Length bare structure	198,2 cm
Width bare structure	89,7 cm
Length structure + H&F Boards + siderails	209,5 cm
Width structure + H&F Boards + siderails	100 cm
Length sleeping surface	197,7 cm
Width sleeping surface	85,8 cm
Variable height minimum	20 cm
Variable height maximum	80 cm
Automatic medium stop position	40 cm
Ground Clearance	24 cm
TYPE	
Variable height type	Y-Type
Type of engine	Electric
Type of brake	Centralized brake in foot only
WEIGHING SYSTEM	
Weighed (Yes/No)	No
MATERIAL & COATING	
Tray Structure	White Epoxy Steel
Base Structure	Black Epoxy Steel
Variable height	Black Epoxy Steel
Bedbase / sleeping surface	Compact HPL
Wheels	Thermoplastic rubber, grey non-marking
BASE	
Type of footing	Castors
Wheels diameter	75 mm
INCLINATIONS	
Backrest	74 °
Legrest	15 °
Knee brake section	38 °
Knee brakes	130 °
Translation	7 cm
Trendelenburg	17 °
Reverse trendelenburg	15 °
PROTECTION SIGN	
Bed	IPX4
Engines	IP66
SAFETY CHARGE	
Maximum patient weight	195 kg
SWL	260 kg
CONFIGURATION	
Configuration	Siderails - H&F boards
Number of lifting pole support	2
Number of IV pole support	2
REGULATION	
Medical device class	Class I
CE marking	CE
European regulation (UE 2017/745)	Yes
FIRE CLASSIFICATION	
Motorization	V0
Compact HPL	M1
Laminated	M1
Melaminated panel	M3

OPTIONS

FUNCTION CONDEMNATION BOX WITH REMOTE-CONTROL OUTLET PLUG & PLAY FOR BEDS Y4
Réf. 4920-50005



BACKUP BATTERY
Réf. 4920-50007



ELECTRIC KNEE BREAK WITH ELECTRIC FOOTREST (V20)
Réf. 4920-50010



SINGLE NIGHTLIGHT POSITIONED UNDER THE BED (LED WITH CONNECTOR)
Réf. 4923-50011



EASYTOUCH SATELLITE WITH THE PRE-PROGRAMMING OF THE MAIN FUNCTIONS (WITHOUT REMOTE CONTROL)
Réf. 4953-50007



PATIENT SATELLITE ON FLEXIBLE ARM - IN ADDITION
Réf. 4953-50008



CAREGIVER SUPERVISOR PLACED AT THE FOOT OF THE BED
Réf. 4953-50015



OPTIONS

BILATERAL RELEASE HANDLE OF BACKREST
Réf. 4920-50008



BILATERAL RELEASE HANDLE OF BACKREST
Réf. 4920-50009



OPTIONS

ALTERNATIVE OPTION ON PACK - MOBILITY OF THE BED TIVOLI 4 EASYMOVE
Réf. 4920-50062



OPTIONS

2 ANNULAR PROTECTIVE BUMPERS FIXED ON THE HEAD & FOOT PANELS
Réf. 4998-00617



ALTERNATIVE OPTION ON PACK - H&F BOARDS VERONE ***LESS VALUE***
Réf. 4998-00712



OPTIONS

NATURAL FINISH OF THE 2 TOP BARS FOR SET OF CLASSIC FOLDABLE SIDE RAILS - 3 BARS
Réf. 4920-50020



WOODEN FINISH OF LONG PANELS - WITHOUT H&F BOARDS
Réf. 4920-50022



UPHOLSTERED FINISH OF LONG PANELS - WITHOUT H&F BOARDS
Réf. 4920-50059



OPTIONS

ALTERNATIVE OPTION ON PACK - WIDTH OF BED TIVOLI 4 120 CM - EMOVE BED LOW
Réf. 4920-50061



OPTIONS

4 RAILS FOR ACCESSORIES FOR TOP OF BED FRAME (2 ON BOTH SIDES)
Réf. 4953-50017



ACCESSORIES

FLEXIBLE HOLDER FOR WIRED REMOTE-CONTROL
Réf. 4920-50049



ACCESSORIES

OPTION HANDRAIL AT THE HEAD FOR SET OF H&F BOARDS VERONE, AUBANCE AND CHEVERNY
Réf. 4998-00713



ALTERNATIVE OPTION ON PACK - SET OF ONE-PIECE H&F BOARDS AUBANCE EASYCLIP Y4 - MELAMINATED DECORS
Réf. 4998-00715



ALTERNATIVE OPTION ON PACK - SET OF ONE-PIECE H&F BOARDS AUBANCE EASYCLIP WITH LACQUERED DECOR Y4 - MELAMINATED DECORS
Réf. 4998-00716



ALTERNATIVE OPTION ON PACK - SET OF WOODEN H&F BOARDS CHEVERNY EASYCLIP - MELAMINATED DECORS
Réf. 4998-00718



ALTERNATIVE OPTION ON PACK - SET OF WOODEN H&F BOARDS CHEVERNY EASYCLIP WITH LACQUERED DECOR - MELAMINATED DECORS
Réf. 4998-00719



ALTERNATIVE OPTION ON PACK - SET OF UPHOLSTERED H&F BOARDS ALPAGA EASYCLIP
Réf. 4998-00721



ALTERNATIVE OPTION ON PACK - H&F BOARDS ATHENA - TWO POSITIONS BOARD AND TRAY
Réf. 4998-00723



ACCESSORIES

SET OF CLASSIC FOLDABLE SIDE RAILS - 3 BARS FOR BED Y4
Réf. 4920-50011



SET OF PADDED PROTECTIONS FOR ¾ SIDE RAILS
Réf. 4923-50023



ALTERNATIVE OPTION ON PACK - SET OF CLASSIC FOLDABLE SIDE RAILS - 3 BARS - FOR BED Y4 ***LESS VALUE***
Réf. ACP-4920-50011



ALTERNATIVE OPTION ON PACK - SET OF CLASSIC FOLDABLE SIDE RAILS WITH ADDITIONAL TOP BAR - 3 BARS (THE PAIR)
Réf. OPP-4920-50012



ALTERNATIVE OPTION ON PACK - SET OF CLASSIC FOLDABLE SIDE RAILS - 4 BARS - ADAPTED TO GROWN CHILDREN (THE PAIR)
Réf. OPP-4920-50015



ACCESSORIES

SET OF WOODEN SIDE RAILS FULL LENGTH - 2 BARS - WITHOUT H&F BOARDS
Réf. 4920-50026



SET OF PADDED PROTECTIONS FOR WOODEN SIDE RAILS FULL LENGTH
Réf. 4920-50060



SET OF ADDITIONNAL TOP BAR FOR WOODEN SIDE RAILS FULL LENGTH
Réf. 4923-50030



ACCESSORIES

SET OF PADDED PROTECTIONS FOR SIDE RAILS ERGOFLEX
Réf. 4491-00017



SET OF ADDITIONNAL TOP BAR FOR ERGOFLEX
Réf. 4491-00018



SET OF ERGOFLEX 60/40 TOTAL CLOSURE - WOODEN FINISH - WITHOUT H&F BOARDS
Réf. 4491-00020



ACCESSORIES

FIXED LIFTING POLE WITH AUTOMATIC REEL
Réf. 4923-50027



LIFTING POLE SWIVELLED BY A PULL WITH REEL - POSITIONED AT THE HEAD OF THE BED
Réf. 4953-50011



4 PLASTIC HOOKS SUPPORT FOR LIFTING POLE
Réf. 4953-50012



ACCESSORIES

SET OF ADDITIONAL BARRIERS FOR CLASSIC SIDE RAILS (THE PAIR)
Réf. 4920-50002



HANDLE TO EASE GETTING OUT OF BED (SOLD PER UNIT) - INCOMPATIBLE WITH KIT XL 110 CM (UNIT)
Réf. 4920-50024



ACCESSORIES

IV POLE WITH FIXED HEIGHT WITH 4 HOOKS
Réf. 4920-50028



IV POLE WITH VARIABLE HEIGHT BY SCREW SYSTEM - 4 HOOKS
Réf. 4953-50013



ACCESSORIES

SUPPORT PLATE FOR REMOVABLE MATTRESS FOR BED LOW EXTENSION
Réf. 4920-50001



ACCESSORIES

VISCOELASTIC FOAM MATTRESS + COVER - WIDTH 120 CM
Réf. 4000-00088



MEMORY FOAM MATTRESS + VISCOELASTIC + COVER FOR KIT XL COMFORT - 200 X 110 X 14 CM - SLOTS ARE CASTED UNDER THE MATTRESS WITHOUT
Réf. 4000-00089



KIT XL COMFORT - SLEEPING PLATFORM 110 CM - INCLUDING : SET OF ONE-PIECE H&F BOARDS AUBANCE WITHOUT LACQUERED STRIP + SLEEPING
Réf. 4920-50046



ACCESSORIES

COMFORTGUARD (UNIT)
Réf. 4491-00019



ACCESSORIES

MATTRESS + COVER - COMFORT - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION
Réf. 4000-00036



MATTRESS + COVER - COMFORT - SLEEPING 20 X 120 CM - FOR BED FRAME EXTENSION
Réf. 4000-00053



MATTRESS + COVER - COMFORT MULTI-ARTICULATION - SLEEPING 200 X 90 CM - RIGHT ANGLES
Réf. 4000-00077



MATTRESS + COVER - COMFORT MULTI-ARTICULATION - SLEEPING 200 X 120 CM - RIGHT ANGLES
Réf. 4000-00080



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 200 X 90 CM - RIGHT ANGLES
Réf. 4000-00081



MATTRESS + COVER - VISCOELASTIC - SLEEPING 200 X 90 CM - RIGHT ANGLES
Réf. 4000-00085



MATTRESS + COVER - NON-MOTORIZED AIR - SLEEPING 200 X 90 CM
Réf. 4000-00092



COVER FOR MATTRESS - WIDTH 90 CM - RIGHT ANGLES
Réf. 4000-00093



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION
Réf. 4000-00098



EXTENSION FOR WAFFLE FOAM MATTRESS + COVER - 20 X 90 X 14 CM
Réf. 4000-00099



MATTRESS + COVER - VISCOELASTIQUE - SLEEPING 20 X 90 CM - FOR BED FRAME EXTENSION
Réf. 4000-00100



MATTRESS + COVER - COMFORT MULTI-ARTICULATION - SLEEPING 200 X 90 CM - TRUNCATED ANGLES
Réf. 4000-00107



MATTRESS + COVER - MULTI-JOINTS WAFFLE - SLEEPING 200 X 90 CM - TRUNCATED ANGLES
Réf. 4000-00111



VISCOELASTIC FOAM MATTRESS + COVER - TRUNCATED ANGLES - WIDTH 90 CM
Réf. 4000-00115



COVER FOR MATTRESS - WIDTH 90 CM - TRUNCATED ANGLES
Réf. 4000-00122



ACCESSORIES

SHOCK ABSORPTION MAT - PROTECTION IN
CASE OF A FALL
Réf. 4000-00065



PATIENT PRESENCE AND/OR POSITION SENSOR
IN THE BED
Réf. 4000-00127



ACCESSORIES

BUMPER ON TOP BED FRAME Y4
Réf. 4920-50031



REMOVABLE AND RETRACTABLE STOPPER FOR
WALL PROTECTION
Réf. 4923-50015



ACCESSORIES

URINAL HOLDER ON SIDE RAILS
Réf. 4923-50013



URINAL POCKET HOLDER
Réf. 4953-50024

